



www.labomed.com
spectro@labomed.com

LB-112 Monocular Biological Microscope

Introduction

With sharp image, competitive and reasonable unit price, LB-112 monocular biological microscopes are ideal instruments for student use. These microscopes are mainly designed for elementary school.

Specifications

Eyepiece:	Huygens Eyepiece H5x, H10x, H16x
Objective:	Achromatic Objective 10x, 40x, 100x
Nosepiece:	Triple Nosepiece
Stage:	Attachable Mechanical Stage 110-120mm / 60x30mm
Condenser:	Abbe with Iris Diaphragm and Filter
Illumination System:	Plan-concave Mirror
Filter:	Blue Filter, Green Filter

