



www.labomed.com
spectro@labomed.com

Introduction

LB-114 Biological Monocular Microscopes are more comfortable and reliable operation with modern design. These microscopes are economical, with all basic features for educational use in elementary and high school.

Applications

These microscopes are economical microscope with all basic features for educational use in elementary and high school.

Technical Specifications

Viewing Head:	Monocular, Inclined at 45°, 360° Rotatable
Nosepiece:	Triple Nosepiece
Eyepiece:	Wide Field Eyepiece WF10×/18
Achromatic Objective:	Achromatic Objective 4×, 10×, 40×(S)
Stage:	Plain Stage with Slide Clips 110×120mm
Focusing:	Coaxial Coarse and Fine Adjustment, Moving Range 10mm
Condenser:	Single Lens NA 0.65 with Iris Diaphragm
Illumination:	LED Illumination, Brightness Adjustable
Package:	Styrofoam & Carton, Dimension 28×20×40cm, 3kg

