



www.labomed.com  
spectro@labomed.com

# LB-132 Portable Measuring Microscope

## Introduction

LB-132 Portable Measuring Microscopes feature smart size, light weight, nice design and easy operation. They are mainly used for measuring or viewing purposes in production or laboratory. They are particularly suitable for spot examination of machinery industry, paper making, printing and textile, etc. They are also suitable to be used by students for education.

## Applications

LB-132 Portable Measuring Microscopes are mainly used for measuring or viewing purposes in production and laboratory area. They are particularly suitable for examination of machinery industry, paper making, printing and textile, etc. They are also suitable to be used by students for education.

## Technical Specifications

Total Magnification:	50×
Mechanical Tube Length:	160mm
Eyepiece:	10×
Objective:	5×
Numerical Aperture:	0.12
Diameter of Viewing field	
in Eyepiece:	φ18mm
Eyepiece Reticle:	0.1mm
Fine Division:	0.02mm
Working Distance:	17.6mm
Focus Range:	30mm
Diameter of Bearing Plane:	φ63mm
Height of Microscope:	210mm
Mass of Microscope:	650g
Illumination:	Pen light, No.7 battery

## Accessories

