





LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-342 Digital LCD Stereo Microscope (5.0MP) with 9 inch HD LCD screen with HD video and HDMI output

9 Inches HD screen (1280×800) real-time display and no need to connect to computer.
Built in 5 Mega digital professional camera and integration of capture, video and storage.
720P/30fp HD video and HDMI Output.
Improving working efficiency and reducing labor fatigue.
It can be connected to a projector and a big monitor screen at the same time.
Images and videos can be stored in the SD card.