




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-401 Trinocular Industrial Inspection and Integrated Circuit Inspection Microscope with Extra Wide Field and Infinity Optical System

Big Stage with Extra Large Moving Range 250×250 mm.

Well designed structure, easy to operate. It is specially used in IT Industry for Large Area Integrated Circuits Circuit, Wafer Observation and Inspection.

Excellent Optical Quality with Infinite Optical System

Stage Speedy Movement by Hand to Get Observation Point Immediately.

Comfortable and Enjoyable Operation with Low and Pre Position Control Knob.

Suitable for different Specimen Observation with Polarization Set.

Easy, Safe and Reliable Lamp Replacement Device.

Available for Measuring with Digital Camera and Software.