



www.labomed.com
spectro@labomed.com

LB-981 Multi-Head Microscopes (for 3 users) with Infinite Optical System (Anti-Fungus)

Introduction

LB-981 Multi-Head Triple View Microscopes (for 3 Users) (Anti-Fungus) are high level microscopes equipped with multi-head for more persons observing at the same time. The optical system is high quality and reliable.

Features

Sharp image display with excellent optical system.

Small space occupancy with integral stand design, low environment requirement with anti-mold technology.

Users friendly and comfortable operation with ergonomic design.

With the features of Infinite Optical System, Effective Illumination, LED Pointer and Image Coherence, it is widely used in Clinic, Scientific Research and Teaching Demonstration.

Applications

LB-981 Multi-Head Triple View Microscopes (for 3 Users) (Anti-Fungus) are mainly used in medical teaching and biological teaching. They also can be used for biological analysis for more experts at the same time. They can be widely used in medical and sanitary establishments, laboratories, institutes, colleges and universities.

Technical Specifications

Optical System:	Infinite Optical System	
Viewing Head:	Seidentopf Trinocular Head, Inclined at 30°, 360° Rotatable, Interpupillary Distance 48-75mm	1Pc
	Seidentopf Binocular Head, Inclined at 30°, 360° Rotatable, Interpupillary Distance 48-75mm	2Pcs
Eyepiece:	Extra Wide Field Eyepiece EW10×/ 20	6Pcs
Nosepiece:	Inverse Quintuple Nosepiece	
Objective:	Infinite Plan Achromatic Objective 4×、10×、40×、100×	
Condenser:	Swing Condenser N.A.0.9/ 0.25	
Focusing:	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm	
Stage:	Double Layers Mechanical Stage 185×142mm, Moving Range 75×55mm	
Kohler Illumination:	External Illumination, Halogen Lamp 12V/100W	
Pointer:	Green LED Pointer, Brightness Adjustable Dual Color Pointer, Brightness Adjustable (optional)	

